

CONTINUOUS PASSIVE MOTION unit is designed for lower limb. It is based on solid-state circuitry, silent, smooth & powerful servo motor for continuous operation, instant start/stop.

APPLICATION

Continuous Passive Motion Therapy if used immediately after surgery can give several important benefits to the patient like maintenance or improvement in the range of motion, reduction in pain and discomfort.



SPECIFICATION: -

Dimension (L*H*D)	1113mm, 240mm, 210mm
Weight	15.5 Kg (Approx.)
Operating Voltage	220v 50Hz
Absorption	60-Watts Max
Fuses	1 Amp
Room Temperature	10° to 40°C
Moisture	10% to 80%
Notation Range	10° to 120°
Digital Timer	01 ~ 99 Minutes Adjustable
Patient Safety	Available