


ADVANCE

AYU-44P-DX

ELECTRO SURGICAL UNIT (400W)



Light Weight

•
Elegant Design

•
Simple Operation

•
The Reliable Economical

•
Digital Display

•
100 Programs

AYU-400W-DX
Electro Surgical

ADVANCE

AYU-44P-DX

ELECTRO SURGICAL UNIT (400W)

High Frequency Electro Surgical Unit with Separate Bipolar Foot Switch

Microprocessor control, innovative design, together with high reliability makes 400W High frequency Surgical Units the most favorable choice for doctors.

In the used surgical diathermy, the technical ability of the unit combined with plain safety are of prime importance –both for the patient and the user, There must be no possibility of accidental , undesirable results in the use of life care systems Electro Surgical Units been developed .

The Ayu 400W-DX is newly developed electro surgical units of particular safety and reliability. The Ayu 400W as the most refined unit of its kind in the world.

Safety in use

The Technical conceptions of the Ayu 400W –DX meet the requirements of highly refined, Electro Surgical techniques all national and international safety standards are conformed to or exceeded. It is not possible for the unit to deliver a higher output than has been set on the control. The micro function allows for an instant reduction in the output for more delicate application such as neurological or micro surgery.

OUTPUT CHARACTERISTICS

	MODE	RATED LOAD (OHMS)	MAXIMUM POWER (WATTS)
Monopolar Cut	PURE	300	400
	BLEND-I	300	300
	BLEND-II	300	250
	ENDO CUT	300	200
Monopolar Coag.	DESICCATE	300	150
	FULGURATE	300	250
	FORCE	300	200
	SPRAY	300	150
	Bipolar	CUT	100
COAG		100	70
MACRO		100	70
MICRO		100	50

Physical Characteristics – LxWxH - 430x405x280 mm
Net Weight - 5Kg. 500Grm. Gross Weight -10Kg. 900Grm.

Standard Accessory



K-708, Sector-23, Sanjay Nagar, Ghaziabad-201002 Delhi NCR, INDIA

Tel : 0120-2784048, 2788618 Cell : +91-9312247417, +91-9654875195

E-mail : info@lcsindia.com, sales@lcsindia.com, anil.singhal22@gmail.com

Website : www.lcsindia.com, www.physiotherapyequipments.com