


ELECTROTHERAPY EQUIPMENTS

LCS 231



CONTINUOUS PASSIVE MOTION OF LCS 231

LCS CPM is designed for lower limb. It is based on solid-state circuitry, silent, smooth & powerful servo motor for continuous operation, instant start/stop.

APPLICATION LCS 231

Continuous Passive Motion therapy immediately used after surgery can give several important benefits to patient like maintenance or improvement in the range of motion, reduction in pain and discomfort.

TRACTION TABLE OF LCS 302

Fixed Height 4 section Traction Table for horizontal Cervical & Lumbar Traction Treatment, comes up with adjustable stool. We offer smooth and comfortable traction for patients.

LCS 302



SPECIFICATION	LCS 231	SPECIFICATION LCS 302
Dimension (L*H*D)	1100mm, 224mm, 185mm	Height 27 Inches (68.58 cm)
Weight	16 Kg. (Approx)	Width 24 Inches (60.96 cm)
Operating Voltage	220V AC 50Hz	Length 94 Inches (238.76 cm)
Absorption	60-Watts Max.	
Fuses	1 Amp	
Room Temperature	10° to 40°C	
Moisture	10% to 80%	
Notation Range	10° to 120°	
Digital Timer	01 ~ 99 Minutes Adjustable	
Patient Safety	Available	

FEATURES LCS 302

- Friction-free lumbar section.
- Comes up with adjustable stool.
- Cushion with comfortable slow recovery foam.
- Lever system to allow front OR back seat angle adjustment.
- LCS 302 table includes a three-piece head section.